





















Icon	Function	Shortcut	Icon	Function	Shortcut
	Axial view	Ctrl + 1		New window	Ctrl + N
	Sagittal view	Ctrl + 2		Close tab	Ctrl + W
	Coronal view	Ctrl + 3		Exit	Ctrl + Q
	Capture current image	Ctrl + S		Local database studies	Ctrl + L
	Capture all of current series	Ctrl + A		PACS	Ctrl + P
	Erase all annotations	Ctrl + E		Open DICOMDIR	Ctrl + D
	Restore viewer to the initial state	Ctrl + R		Open files	Ctrl + O
	Cobb angle	⇧ + A		Open files from a directory	Ctrl + ⇧ + D
	Polyline ROI	⇧ + R		Maximise to multiple screens	Ctrl + ⇧ + M
	Magnify	⇧ + Z		Move to previous screen	Ctrl + ⇧ + ←
	Desynchronise all viewers	⇧ + U		Move to next screen	Ctrl + ⇧ + →
	External applications	⇧ + F1 ... F12			